

NFAT2 Rabbit pAb


ABZ1443

Product Overview


Name	NFAT2 Rabbit pAb
Accession(Primary)	O95644
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 101 kDa; Observed MW: 105 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:500; IF/ICC1:100-1:500; ELISA-1:5000-1:20000

Target information (O95644)

Gene Name	NFATC1
About Gene Symbol	



Western blot analysis of NFAT2 in HuvEc lysates using NFAT2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung using NFAT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.