


GABARAP Rabbit pAb

ABZ1444


Product Overview

Name	GABARAP Rabbit pAb
Accession(Primary)	O95166
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14,16 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (O95166)

Gene Name GABARAP**About Gene Symbol**

Western blot analysis of GABARAP in A549, 3T3 and HeLa lysates using GABARAP antibody.



Western blot analysis of GABARAP in rat Hippocampus, rat Cerebral Cortex, rat Brain and A549 lysates using GABARAP antibody.