

CD133 Rabbit pAb


ABZ1446

Product Overview

Name	CD133 Rabbit pAb
Accession(Primary)	O43490
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 97 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000

Target information (O43490)

Gene Name	PROM1
About Gene Symbol	



Western blot analysis of CD133 in A549, SW620, HT-29, Jurkat, mouse kidney, mouse lung lysates using CD133 antibody.