

FRA1 Rabbit pAb


ABZ1449

Product Overview


Name	FRA1 Rabbit pAb
Accession(Primary)	P15407
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P15407)

Gene Name	FOSL1
About Gene Symbol	



Western blot analysis of FRA1 in HeLa, Jurkat
lysates using FRA1 antibody



Immunohistochemistry analysis of paraffin-
embedded Human gastric cancer using FRA1
antibody. High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.