


## UNC5CL Rabbit pAb

ABZ1450

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | UNC5CL Rabbit pAb                          |
| <b>Accession(Primary)</b>   | Q8IV45                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 58 kDa; Observed MW: 58 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000                            |

### Target information (Q8IV45)

**Gene Name** UNC5CL**About Gene Symbol**

Western blot analysis of UNC5CL in mouse kidney,  
293T lysates using UNC5CL antibody

