

GPC4 Rabbit pAb

ABZ1451

Product Overview

Name	GPC4 Rabbit pAb
Accession(Primary)	O75487
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 62 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:2000


Target information (O75487)

Gene Name GPC4**About Gene Symbol**

kDa

130—
95—
72—
55—

36—
28—
17—



Western blot analysis of GPC4 in Human pancreas
lysates using GPC4 antibody

