

HIF1 alpha Rabbit pAb

ABZ1452

Product Overview


Name	HIF1 alpha Rabbit pAb
Accession(Primary)	Q16665
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 92-130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; IP-1:200; ELISA-1:40000

Target information (Q16665)


Gene Name	HIF1A
About Gene Symbol	



Western blot analysis of HIF1 alpha in LOVO lysates using HIF1 alpha antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon tissue using HIF1 alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat lung tissue using HIF1 alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.