

FSHR Rabbit pAb


ABZ1454

Product Overview

Name	FSHR Rabbit pAb
Accession(Primary)	P23945
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; ELISA-1:5000; IF-1:50-1:200

Target information (P23945)

Gene Name	FSHR
About Gene Symbol	



Western blot analysis of FSHR in HeLa lysates
using FSHR antibody.