

TSG101 Rabbit pAb


ABZ1457

Product Overview


Name	TSG101 Rabbit pAb
Accession(Primary)	Q99816
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:5000

Target information (Q99816)

Gene Name	TSG101
About Gene Symbol	



Western blot analysis of TSG101 in COS7 lysates using TS101 antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of TSG101 in HeLa using TGS101 antibody (green)