

ZAP70 Rabbit pAb


ABZ1458

Product Overview

Name	ZAP70 Rabbit pAb
Accession(Primary)	P43403
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P43403)

Gene Name	ZAP70
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue, using ZAP70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.