

ADCY5/6 Rabbit pAb


ABZ1459

Product Overview


Name	ADCY5/6 Rabbit pAb
Accession(Primary)	O95622
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 139 kDa; Observed MW: 139 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (O95622)


Gene Name	ADCY5/ADCY6
About Gene Symbol	




Western blot analysis of ADCY5/6 in various lysates using ADCY5/6 antibody.



Western blot analysis of ADCY5/6 in COLO205 lysates using A Cyclase V/VI antibody.



Western blot analysis of ADCY5/6 in COLO25 lysates using ADCY5/6 antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of ADCY5/6 in HUVEC lysates using ADCY5/6 antibody.