

CNGB1 Rabbit pAb

ABZ1462


Product Overview


Name	CNGB1 Rabbit pAb
Accession(Primary)	Q14028
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 140 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q14028)


Gene Name CNGB1

About Gene Symbol





Western blot analysis of CNGB1 in 3T3 lysates using CNGB1 β antibody.



Western blot analysis of CNGB1 in NIH/3T3 and A549 lysates using CNGB1 antibody. The lane on the right is blocked with the synthesized peptide.