

EGR2 Rabbit pAb


ABZ1463

Product Overview


Name	EGR2 Rabbit pAb
Accession(Primary)	P18146
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P18146)

Gene Name	EGR2
About Gene Symbol	



Western blot analysis of EGR2 in various lysates using Egr2 antibody.



Western blot analysis of EGR2 in HeLa and Jurkat lysates using EGR1/2 antibody. The lane on the right is blocked with the synthesized peptide.