

HEXIM1 Rabbit pAb


ABZ1465

Product Overview


Name	HEXIM1 Rabbit pAb
Accession(Primary)	O94992
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (O94992)


Gene Name	HEXIM1
About Gene Symbol	HEXIM1 (HEXIM P-TEFb Complex Subunit 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cardiovascular disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HEXIM1 in rat HEART lysates using HEXIM1 antibody.



Immunofluorescence analysis of HEXIM1 in HeLa using HEXIM1 antibody. The picture on the right is blocked using the synthesized peptide.



Western blot analysis of HEXIM1 in NIH/3T3 lysates using HEXIM1 antibody. The lane on the right is blocked with the synthesized peptide.

