

## Laminin alpha 5 Rabbit pAb


ABZ1466

### Product Overview


|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | Laminin alpha 5 Rabbit pAb                                          |
| <b>Accession(Primary)</b>   | O15230                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse                                                        |
| <b>Calculated MW</b>        | Calculated MW: 400 kDa; Observed MW: 400 kDa                        |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC:-1:100-1:300;<br>IF-1:100-1:300; ELISA-1:20000 |

### Target information (O15230)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | LAMA5 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of Laminin alpha 5 in RAW, HepG2, A549 lysates using Laminin alpha 5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung carcinoma tissue, using LAMA5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.