

MYOM2 Rabbit pAb

ABZ1469


Product Overview

Name	MYOM2 Rabbit pAb
Accession(Primary)	P54296
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 165 kDa; Observed MW: 165 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300


Target information (P54296)

Gene Name MYOM2


About Gene Symbol



Western blot analysis of Myomesin 2 in various lysates using Myomesin 2 antibody.



Western blot analysis of Myomesin 2 in HT-29 lysates using Myomesin2 antibody.



Western blot analysis of Myomesin 2 in COS7 and HT-29 lysates using MYOM2 antibody. The lane on the right is blocked with the synthesized peptide.