

PDPK1 Rabbit pAb


ABZ1470

Product Overview


Name	PDPK1 Rabbit pAb
Accession(Primary)	O15530
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (O15530)

Gene Name	PDPK1
About Gene Symbol	PDPK1 (3-Phosphoinositide Dependent Protein Kinase 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely stu



Western blot analysis of PDK1 in COLO205 lysates using PDK1 antibody.



Western blot analysis of PDK1 in HELA COLO 293T lysates using PDK1 antibody.