

PKNOX2 Rabbit pAb


ABZ1471

Product Overview

Name	PKNOX2 Rabbit pAb
Accession(Primary)	Q96KN3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 52 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q96KN3)

Gene Name	PKNOX2
About Gene Symbol	



Western blot analysis of PKNOX2 in HepG2 lysates using PKNOX2 antibody.