

Tau Rabbit pAb


ABZ1472

Product Overview


Name	Tau Rabbit pAb
Accession(Primary)	P10636
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P10636)


Gene Name	MAPT
About Gene Symbol	



Western blot analysis of Tau in various lysates using Tau antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Tau antibody. High-pressure and temperature Tris-EDTA pH 8.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Tau in mouse brain using Tau antibody. The lane on the right is blocked with the synthesized peptide.