

TROP2 Rabbit pAb


ABZ1473

Product Overview


Name	TROP2 Rabbit pAb
Accession(Primary)	P09758
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P09758)

Gene Name	TACSTD2
About Gene Symbol	TACSTD2 (Tumor Associated Calcium Signal Transducer 2) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify the



Western blot analysis of TACSTD2 in various lysates using TACSTD2 antibody.



Western blot analysis of TACSTD2 in A549 and NIH/3T3 lysates using TACD2 antibody. The lane on the right is blocked with the synthesized peptide.