

ZFP106 Rabbit pAb

ABZ1474


Product Overview

Name	ZFP106 Rabbit pAb
Accession(Primary)	Q9H2Y7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 209 kDa; Observed MW: 209 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000


Target information (Q9H2Y7)

Gene Name ZNF106

About Gene Symbol



Western blot analysis of ZFP106 in COS7 lysates using ZFP106 antibody.



Western blot analysis of ZFP106 in COS7 lysates using ZFP106 antibody.