

Histone H1.8 Rabbit pAb

ABZ1476


Product Overview

Name	Histone H1.8 Rabbit pAb
Accession(Primary)	P16401
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (P16401)

Gene Name H1-8

About Gene Symbol



Immunofluorescence analysis of Histone H1.8 in HeLa cells using Histone H1.8 antibody. Sample with blocking peptide on the right.