



Cytochrome P450 4Z1 Rabbit pAb

ABZ1478


Product Overview

Name	Cytochrome P450 4Z1 Rabbit pAb
Accession(Primary)	Q86W10
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (Q86W10)

Gene Name CYP4Z1**About Gene Symbol**

Immunofluorescence analysis of Cytochrome P450 4Z1 in HUVEC cells using Cytochrome P450 4Z1 antibody. Sample with blocking peptide on the right.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Cytochrome P450 4Z1 antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.