


SGLT1 Rabbit pAb

ABZ1479

Product Overview

Name	SGLT1 Rabbit pAb
Accession(Primary)	P13866
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P13866)

Gene Name SLC5A1**About Gene Symbol**

Western blot analysis of SGLT1 in HepG2 lysates using SGLT1 antibody.

