

AMBN Rabbit pAb

ABZ1480

Product Overview

Name	AMBN Rabbit pAb
Accession(Primary)	Q9NP70
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q9NP70)

Gene Name AMBN**About Gene Symbol**

