

NCAPG2 Rabbit pAb

ABZ1481


Product Overview

Name	NCAPG2 Rabbit pAb
Accession(Primary)	Q86XI2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q86XI2)

Gene Name NCAPG2

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human lung carcinoma, using NCAPG2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.