

Rab2B Rabbit pAb

ABZ1482


Product Overview


Name	Rab2B Rabbit pAb
Accession(Primary)	Q8WUD1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q8WUD1)

Gene Name RAB2B

About Gene Symbol





Western blot analysis of Rab2B in 293 lysates
using RAB2B antibody.