


Ubinuclein 1 Rabbit pAb

ABZ1483


Product Overview

Name	Ubinuclein 1 Rabbit pAb
Accession(Primary)	Q9NPG3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (Q9NPG3)

Gene Name UBN1**About Gene Symbol**

Immunofluorescence analysis of Ubinuclein 1 in A549 cells using Ubinuclein 1 antibody. Sample with blocking peptide on the right.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Ubiquitin 1 antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.