

GAS41 Rabbit pAb

ABZ1484

Product Overview

Name	GAS41 Rabbit pAb
Accession(Primary)	O95619
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (O95619)

Gene Name YEATS4**About Gene Symbol****(kD)**

117-

85-

48-

34-

26-

19-



GAS41

Western blot analysis of GAS41 in HT-29 lysates
using GAS41 antibody.

