

PEX16 Rabbit pAb

ABZ1485

Product Overview

Name	PEX16 Rabbit pAb
Accession(Primary)	Q9Y5Y5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 39 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000-1:20000

Target information (Q9Y5Y5)

Gene Name PEX16**About Gene Symbol**

