


NT5DC3 Rabbit pAb

ABZ1486

Product Overview

Name	NT5DC3 Rabbit pAb
Accession(Primary)	Q86UY8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 63 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000

Target information (Q86UY8)

Gene Name NT5DC3**About Gene Symbol**

Western blot analysis of NT5DC3 in A431 lysates using NT5DC3 antibody.

