

ZBTB45 Rabbit pAb

ABZ1489


Product Overview

Name	ZBTB45 Rabbit pAb
Accession(Primary)	Q96K62
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500; IF-1:50-1:200


Target information (Q96K62)

Gene Name ZBTB45

About Gene Symbol



Western blot analysis of ZNF499 in Hela, Jurkat, mouse Brain, mouse Kidney, rat Brain lysates using ZNF499 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse Kidney Tissue using ZNF499 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.