

Lactate Dehydrogenase D Rabbit pAb

ABZ1490


Product Overview

Name	Lactate Dehydrogenase D Rabbit pAb
Accession(Primary)	Q86WU2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q86WU2)

Gene Name LDHD

About Gene Symbol



Western blot analysis of Lactate Dehydrogenase D in Hela, mouse Skeletal Muscle, mouse Liver, rat Kidney lysates using LDHD antibody.