

Pepsinogen II/PGC Rabbit pAb


ABZ1491

Product Overview

Name	Pepsinogen II/PGC Rabbit pAb
Accession(Primary)	P20142
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Plants
Calculated MW	Calculated MW: 42 kDa; Observed MW: 105-109 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (P20142)

Gene Name	PGC
About Gene Symbol	



Western blot analysis of PEPC in Maize lysates
using PEPC antibody.