

## SLUG Rabbit pAb

ABZ1492


### Product Overview

|                             |                                                 |
|-----------------------------|-------------------------------------------------|
| <b>Name</b>                 | SLUG Rabbit pAb                                 |
| <b>Accession(Primary)</b>   | O43623                                          |
| <b>Host Species</b>         | Rabbit                                          |
| <b>Clonality</b>            | Polyclonal                                      |
| <b>Isotype</b>              | IgG                                             |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                               |
| <b>Calculated MW</b>        | Calculated MW: 30 kDa; Observed MW: 30 kDa      |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:200-1:500; IF-1:50-1:200 |


### Target information (O43623)

**Gene Name** SNAI2

**About Gene Symbol**



Western blot analysis of SNAI2 in MCF-7, mouse Heart, rat Heart, rat Brain lysates using SNAI2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Breast carcinoma using SNAI2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.