


TLR4Rabbit pAb

ABZ1493

Product Overview

Name	TLR4Rabbit pAb
Accession(Primary)	O00206
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:200-1:500; IF-1:50-1:200

Target information (O00206)

Gene Name TLR4**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded mouse Spleen Tissue using Toll-Like Receptor 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.