

CACNG3 Rabbit pAb

ABZ1495

Product Overview

Name	CACNG3 Rabbit pAb
Accession(Primary)	O60359
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC-1:200-1:500; IF-1:50-1:200

Target information (O60359)

Gene Name	CACNG3
About Gene Symbol	



Western blot analysis of CACNG3 in Human Brain, mouse Brain, rat Brain lysates using CACNG3 antibody.