

IGFN1 Rabbit pAb


ABZ1496

Product Overview

Name	IGFN1 Rabbit pAb
Accession(Primary)	Q86VF2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 138 kDa; Observed MW: 115 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-1:200-1:500; IF-1:50-1:200

Target information (Q86VF2)

Gene Name	IGFN1
About Gene Symbol	



Western blot analysis of Elongation Factor 1A2 Binding Protein in Pig Skeletal Muscle, rat Skeletal Muscle, rat Skeletal Muscle lysates using eEF1A2 antibody.