

Cav gamma 4 Rabbit pAb


ABZ1497

Product Overview

Name	Cav gamma 4 Rabbit pAb
Accession(Primary)	Q9UBN1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 37 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC-1:200-1:500; IF-1:50-1:200

Target information (Q9UBN1)

Gene Name	CACNG4
About Gene Symbol	



Western blot analysis of Cav gamma 4 in Human Brain, mouse Brain, rat Brain lysates using Cav gamma 4 antibody.