

Cav gamma 5 Rabbit pAb

ABZ1498


Product Overview

Name	Cav gamma 5 Rabbit pAb
Accession(Primary)	Q9UF02
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200; IF-1:50-1:200

Target information (Q9UF02)

Gene Name CACNG5

About Gene Symbol



Western blot analysis of Cav gamma 5 in Human Brain, mouse Brain, rat Brain lysates using Cav gamma 5 antibody.