

Superoxide Dismutase 1 Rabbit pAb

ABZ1499


Product Overview

Name	Superoxide Dismutase 1 Rabbit pAb
Accession(Primary)	P00441
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:1:50-1:300; IF-1:50-1:200

Target information (P00441)

Gene Name SOD1

About Gene Symbol



Western blot analysis of Superoxide Dismutase 1 in HeLa, MCF-7, mouse Brain, rat Brain lysates using SOD1 antibody.