

GAL4 Activation Domain Rabbit pAb

ABZ1500


Product Overview

Name	GAL4 Activation Domain Rabbit pAb
Accession(Primary)	P04386
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Saccharomyces cerevisiae
Calculated MW	Refer to figures
Recommended Dilution	WB-1:10000-1:50000

Target information (P04386)

Gene Name GAL4

About Gene Symbol



Western blot analysis of GAL4 Activation Domain in Recombinant AD protein lysates using GAL4 Activation Domain antibody.

