

LexA DNA Binding Region Rabbit pAb

ABZ1502


Product Overview

Name	LexA DNA Binding Region Rabbit pAb
Accession(Primary)	P0A7C2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Cyanobacteria
Calculated MW	Refer to figures
Recommended Dilution	WB-1:10000-1:50000

Target information (P0A7C2)

Gene Name lexA

About Gene Symbol



Western blot analysis of LexA DNA Binding Region in Recombinant LexA protein lysates using LexA DNA Binding Region antibody.

