

GPR75 Rabbit pAb

ABZ1503


Product Overview


Name	GPR75 Rabbit pAb
Accession(Primary)	O95800
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:200-1:1000

Target information (O95800)

Gene Name GPR75

About Gene Symbol





Immunofluorescence analysis of GPR75 in LOVO using GPR75 antibody. The picture on the right is blocked using the synthesized peptide.