

NOX2/gp91phox Rabbit pAb

ABZ1504

Product Overview


Name	NOX2/gp91phox Rabbit pAb
Accession(Primary)	P04839
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 65 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P04839)


Gene Name	CYBB
About Gene Symbol	




Western blot analysis of NOX2 in K562 lysates using NOX2 antibody.




Western blot analysis of NOX2 in HEPG2 lysates using NOX2 antibody.




Immunohistochemistry analysis of paraffin-embedded Human lymph node using NOX2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lymph node using NOX2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lymph node using NOX2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human liver using NOX2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.