


Apelin Rabbit pAb

ABZ1516


Product Overview

Name	Apelin Rabbit pAb
Accession(Primary)	Q9ULZ1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:100-1:300; ELISA-1:20000


Target information (Q9ULZ1)

Gene Name APLN**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded rat brain using Apelin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Apelin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat brain using Apelin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.