

## NTE Rabbit pAb

ABZ1521

### Product Overview


|                             |                                                                    |
|-----------------------------|--------------------------------------------------------------------|
| <b>Name</b>                 | NTE Rabbit pAb                                                     |
| <b>Accession(Primary)</b>   | Q8IY17                                                             |
| <b>Host Species</b>         | Rabbit                                                             |
| <b>Clonality</b>            | Polyclonal                                                         |
| <b>Isotype</b>              | IgG                                                                |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                  |
| <b>Calculated MW</b>        | Calculated MW: 151 kDa; Observed MW: 151 kDa                       |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC:-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:20000 |

### Target information (Q8IY17)


|                          |        |
|--------------------------|--------|
| <b>Gene Name</b>         | PNPLA6 |
| <b>About Gene Symbol</b> |        |




Western blot analysis of NTE in AD293 lysates using NTE antibody.



Immunohistochemistry analysis of paraffin-embedded rat testis using NTE antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat testis using NTE antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat testis using NTE antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.