

FMO3 Rabbit pAb

ABZ1524

Product Overview


Name	FMO3 Rabbit pAb
Accession(Primary)	P31513
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P31513)


Gene Name	FMO3
About Gene Symbol	




Western blot analysis of FMO3 in HeLa lysates using FMO3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using FMO3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human liver using FMO3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lung using FMO3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.