

Collagen V alpha 2 Rabbit pAb


ABZ1526

Product Overview


Name	Collagen V alpha 2 Rabbit pAb
Accession(Primary)	P05997
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 145 kDa; Observed MW: 145 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P05997)


Gene Name	COL5A2
About Gene Symbol	




Western blot analysis of Collagen V alpha 2 in HepG2 lysates using Collagen V alpha 2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human pancreas using Collagen V alpha 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using Collagen V alpha 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using Collagen V alpha 2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.