

C4A/B Rabbit pAb


ABZ1527

Product Overview


Name	C4A/B Rabbit pAb
Accession(Primary)	P0C0L4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 193 kDa; Observed MW: 180 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P0C0L4)


Gene Name	C4A/C4B
About Gene Symbol	




Western blot analysis of C4A/B in AD293 lysates using C4A/B antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using C4A/B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human liver using C4A/B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using C4A/B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.